



Dyna-Form™ Trolley Mattress

The **Dyna-Form Trolley Mattress** has been designed especially for users considered to be that of “**High Risk**” of pressure ulcer development in potentially lengthy A&E waiting environments. This product can be made to measure in order to fit most trolley designs. Dyna-Form design quality comprises of a specialist “castellated cut” CME foam surface, complete with a high frequency welded, multi-stretch, vapour permeable cover and non-slip base. This product is ideal for use in busy A&E departments that have an extensive patient turnover.

Features

①	CME pressure reducing foam insert	Provides patient comfort, whilst maintaining excellent pressure re-distribution
②	Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear and friction forces, whilst assisting infection control policies
③	High frequency welded cover	Prevents fluid ingress & aids decontamination processes (can be laundered up to 95°C)
④	Non-slip base	Prevents slippage, maintaining patient safety
⑤	Fully replaceable parts	Ensures great value for money

Benefits

Technical Specifications

Product Code	TROLLEY/MAT/SIZE
Risk Category	High Risk
Weight Limit	26 Stone / 165kg
Warranty	1 Year Foam / 1 Year Cover
Dimensions	Made to measure
Foam Density	34–36kg/m ³ Nominal Hardness – 125 Newtons
Fatigue Class	Severe
Fire Retardancy	BS7177; 2008 (Crib 5 Composite)
Weight	Approximately 1kg
Colour	Blue
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)